



高鐵儀器 精密卓越



GOTECH "the best partner" to raise your competition

GT-7017-URP 換氣式老化箱

Aging Oven

IEC-60811

本機適用於試驗皮類、塑膠類、橡膠類、布類等材料受熱前與受熱後變化特性。根據調節流量調節閥,以流量測定送入試驗槽的空氣流量,來計算出換氣率。本機在規定之溫度及時間下均勻受熱後,觀測試樣老化前後之耐黃、開膠、收縮、伸長、殘餘率等性能,從而判定其老化特性。

GT-7017-URP aging oven is designed to determine the characteristics of leather, plastics, rubber, textile before and after heating. This tester is equipped with flow control valve which allows the operator to adjust the air change rate. Test specimens can be heated uniformly at controlled temperature until the specified period of time is reached.



加熱方式 Heating system	熱風循環(電熱式SSR控制) Air circulation and SSR
控制方式 Control system	HMI+PLC
溫度範圍 Temp. range	常溫 Ambient temp. +10°C ~ 300°C
溫度最小讀數 Resolution	0.1°C
換氣率範圍 Air exchange rate	(3 ~ 50)次/小時 times/hour
回轉盤速度 Turntable rotating speed	(5 ~ 10)r.p.m.
內箱尺寸 Interior dimension(W×D×H)	45×45×50 cm
體積 Exterior dimension(W×D×H)	120×68×165 cm
重量(約) Weight (approx.)	296 kg
電源 Power	1φ,AC 220V (依國別或指定 or specify)

瞭解更多產品資訊 · 請登錄我司網站 www.gotech.biz



創於1977年

高鐵科技股份有限公司

GOTECH TESTING MACHINES INC.

本廠: 台中市工業區24路28號

二廠: 台中市工業區24路31號

TEL: (04)23590061

FAX: (04)23591553

網址: www.gotech.biz E-mail: gt112@gotech.biz

高鐵檢測儀器(東莞)有限公司

GOTECH TESTING MACHINES CO., LTD.

地址: 廣東省東莞市南城區宏圖高科技開發區宏三路

TEL: 0769-89986666 22400899 FAX: 0769-22400699

E-mail: gt-d130@gotech.biz

手機: 13926831580

高特威爾檢測儀器(青島)有限公司

GOODTECHWILL TESTING MACHINES CO., LTD.

地址: 中國青島市市北區哈爾濱路52號青建太陽島1號商務樓8層

TEL: 0532-85931956 85931957

FAX: 0532-85931950

E-mail: gt-q101@gotech.biz

手機: 13708970121