



# 高鐵儀器 精密卓越



## GOTECH "the best partner" to raise your competition

### GT-7008-TR 低溫回縮試驗機

### Low Temperature Retraction Tester

GB/T 7758-2020 · ISO 2921:2019 · ASTM D1329-2016 · JIS K 6261-2006 7.0

本試驗機用於測定拉伸的硫化橡膠經冷卻後的溫度回縮性能。將試樣在標準實驗室溫度下進行拉伸，然後冷卻到在除去拉伸力時不會出現回縮的溫度。除去拉伸力，並以均勻的速率升高溫度。測定出現規定回縮率時的溫度。

TR tester is designed for determining the retraction characteristics of the stretched vulcanized rubber after cooling. The specimen is stretched at a standard laboratory temperature, and then is cooled to a temperature which retraction will not occur after removing tensile force. The tensile force is removed and temperature is increased at a constant rate. The temperature which the specified retraction rate occurs is determined.



溫度範圍 Temp. range	常溫~ -70°C Ambient temp. to -70°C
計數組數 Test station	3 組/6 組【選一】3 or 6 (please specify)
溫度精度 Temp. resolution	0.1°C
溫度偏差 Temp. accuracy	±2.0°C
升溫速率 Heating rate	1.0°C/min
製冷方式 Cooling system	壓縮機製冷 Compressor cooling system
最大伸長率 Maximum extension range	350%(依50mm試樣計算)
位移精度 Displacement resolution	0.01mm
控制方式 Control system	低溫回縮專用電腦軟體Dedicated software
體積 (W×D×H) Dimension (W×D×H)	105×85×235 cm
電源 Power	1φ AC 220V 50Hz (依國別或指定 please specify)

瞭解更多產品資訊 · 請登錄我司網站 [www.gotech.biz](http://www.gotech.biz)



創於1977年

高鐵路科技股份有限公司

GOTECH TESTING MACHINES INC.

本廠: 台中市工業區24路28號

二廠: 台中市工業區24路31號

TEL: (04)23590061

FAX: (04)23591553

網址: [www.gotech.biz](http://www.gotech.biz) E-mail: [gt112@gotech.biz](mailto:gt112@gotech.biz)

高鐵路檢測儀器(東莞)有限公司

GOTECH TESTING MACHINES CO., LTD.

地址: 中國廣東省東莞市南城街道宏三路2號

TEL: 0769-89986666 22400899 FAX: 0769-22400699

E-mail: [gt-d130@gotech.biz](mailto:gt-d130@gotech.biz) 手機: 13926831580

高特威爾檢測儀器(青島)有限公司

GOODTECHWILL TESTING MACHINES CO., LTD.

地址: 中國山東省青島市城陽區招商網穀合心園19號樓101

TEL: 0532-85931956 85931957 FAX: 0532-85931950

E-mail: [gt-q101@gotech.biz](mailto:gt-q101@gotech.biz) 手機: 13708970121